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
012167 RE-1B Reference electrode (Ag/AgCl)

• Detaching the protective cap

(Be careful that a transparent heat shrinkable tube would be separated from the glass tube if you pull the black protective cap in too strong force.)


CAUTION:

1. Do not remove the label of electrode.
The label is required for any case of inquiry.
2. Refrain the electrode from using in strong acid or base solution.
3. Use the electrode at room temperature and atmospheric pressure.
4. To avoid breaking the glass tube, the electrode should be protected from strong shock.
5. The porous glass of tip may be discolored after use, which is due to the nature of the glass. Hence, we can not respond to exchange.
6. The inner solution can not be changed.
7. This electrode is intended for aqueous. Do not use it in organic solvent.

• Attention for setup

If air bubbles were present around the porous glass, slightly flick the electrode to clear the bubbles. The air bubbles may obstruct the liquid conduction between inner solution and external solution, that may finally cause the electrode potential abnormal.


• How to keep the reference electrode

The inner solution of the electrode is 3 M NaCl solution. To stabilize the electrode potential, 3 M NaCl solution should be used at keeping.

If the electrode does not used for long time, it should be stored in our preservative vial (optional) to avoid vaporizing of inner solution. If the electrode is kept in different chloride ion concentration solution, the electrode potential may not be kept correctly due to the chloride ion concentration change of inner solution.

Clean the electrode with distilled water after use to avoid contamination.

012108 RE-PV
Preservative vial for Reference electrode

You can browse the checking data of electrode in below URL.

<http://www.als-japan.com/dl/>

Manufactured by: ALS Co., Ltd
URL: <http://www.als-japan.com>
e-mail: sales@als-japan.com

Exclusive distributor: BAS Inc.